



Density Determination Kit

KDK-300 is suitable for measuring density of solids as well as liquid samples.



MODEL: KDK-300

KDK-300 is suitable for all weighing balance if pan size is less than \varnothing 125 mm. If your weighing balance does not have a density function, user can calculate the density using the below formula.

$$\text{Solid Density} = \frac{(\text{Weight in Air}) \times (\text{SG of fluid})}{(\text{gm/cm}^3) \quad \text{Weight in Air} - \text{Weight in fluid}}$$

For liquid density measurement, first measure volume of crystal ball in distilled water. Then use the volume obtained here when measuring density of test liquid.

$$\text{Volume of Crystal} = \frac{\text{Weight in Air} - \text{Weight in Liquid}}{0.997\text{g/cm}^3(\text{Density of water}) - 0.001185\text{ g/cm}^3(\text{Buoyaoncy of Air})}$$

(cm³)

$$\text{Liquid Density} = \frac{\text{Weight in Air} - \text{Weight in Liquid}}{\text{Volume of Crystal Ball}}$$

(gm/cm³)

STANDARD

ACCESSORIES:

- Beaker
- Stand for beaker
- Hanging Holder
- Object Holder
- Crystal Ball

LABORATORY & INDUSTRIAL
SOLUTIONS