



Vacuum Filtration Manifolds



Vacuum Filtration Manifold

C



Vacuum Filtration Manifold.


Materials:

- Body and branches : Type 316L stainless steel
- Valves : PTFE stopcock, Electroplated SS body.

Connections:

- Filter Head : 6 Stage, 3 Stage, Single Stage.
- Supported Funnel : Glass , PP And Stainless Steel.

Product Code	Description
VFM01	Vacuum Filtration Manifold 1 Stage
VFM03	Vacuum Filtration Manifold 3 Stage
VFM06	Vacuum Filtration Manifold 6 Stage



Vacuum Filtration Manifold GMP MODEL.

Materials:

- Body and branches : Type 316L stainless steel
- Valves : PTFE stopcock, Electroplated SS body.

Connections:

- Filter Head : 6 Stage, 3 Stage, Single Stage.
- Supported Funnel: Disposable Cups

Product Code	Description
GMP01	Vacuum Filtration Manifold 1 Stage
GMP03	Vacuum Filtration Manifold 3 Stage
GMP06	Vacuum Filtration Manifold 6 Stage



Microfilt 47

Materials:

- Body and branches : Type 316L stainless steel
- Valves : PTFE stopcock, Electroplated SS body.

Connections:

- De-attachable Filter Head
- Filter Head : 6 Stage, 3 Stage, 1 Stage
- Supported Funnel : Microfilt funnels.

Product Code	Description
MF47-01	Mircofilt 47 (1 Stage)
MF47-03	Mircofilt 47 (3 Stage)
MF47-06	Mircofilt 47 (6 Stage)



Microfilt 47 Complete Kit

In the Kit:

- 3/6 Stage Manifold
- Micro Filt Disposable filter Funnel
- Oil Free Vacuum Pump AE45s
- 47mm x 0.45um Cellulose Nitrate Membrane Filter Paper
- Vacuum Trap with Silicon Tube
- Forceps

Product Code	Description
MF47C-01	Microfilt 47 Complete Kit (1 Stage)
MF47C-03	Microfilt 47 Complete Kit (3 Stage)
MF47C-06	Microfilt 47 Complete Kit (6 Stage)