



Safety Products & Accessories Z



Hand gloves Dispenser

Made Of ss 304
Light weight & Wallmounted
Available small , Midium & Large Size (Customised Size Available)



Hairnet Dispenser

Made Of Clear Acrylic
Light weight & Wall mounted
Available small , Medium & Large Size (Customized Size Available)



Syringe Filter & HPLC Vial Dispenser

No. Of Housing: 2 , 3 , 4 Station For Store HPLC Vials & Syringe filter storage Easy front rotary Knobs to collect sample untouched
Material: Acrylic Housing & Body SS 304



HPLC Vials hplc screw vials with insert cap and septum

Vial Volume: 1.5-2.0ml
Vial Material: Borosilicate Glass/Glass Type 1
Application: Lab consumables for analysis
Certification: ISO9001, SGS, COA
Vial Thread: Screw
Cap: Plastic open/solid cap
Septa: Non-slit, Pre-Slit, PTFE/Silicone



HPLC Vial Rack

Wells:50 Holes
Material: Made by polypropylene
Design for 1.5ml/2ml HPLC autosampler vial, Aijiren can also supply other vial racks for 4ml, 8ml, 10ml, 20ml, 30ml, 40ml
Easy handing and transportation of sample vial;
vial rack can be stacked in the empty or filled with a set of vials
Stable standing position because of solid construction;
Alphanumeric and indexing are for easy tracking



HPLC Last Drop Filter

Last Drop Inlet Filters have a unique inlet design that allows your system to use almost every "drop" of mobile phase from the
for 1/16th and 1/8th tubing we suggest flange free connectors.