



Vacuum Filtration Unit



Vacuum Filter Holder

A

Glass Filter Holder

- MOC : Borosilicate Glass
- Support screen : Sintered Glass
- Stopper : Silicone Rubber
- Clamp : Anodized Aluminum
- Filter Size : 25mm, 47mm, 90 mm
- Filtration area : 2.1 cm², 9.6 cm², 43 cm²
- Funnel Capacity: 300ml, 500ml, 1000ml.



Product Code	Description
GFH300	Glass Filter Holder 300ml
GFH500	Glass Filter Holder 500ml
GFH1000	Glass Filter Holder 1000ml

Stainless Steel Filter Holder

- MOC : SS316
- Support screen : SS316
- Stopper : Silicone Rubber
- Clamp : Inbuilt SS316
- Filter Size : 47mm,
- Filtration area : 13.6 cm²
- Funnel Capacity: 300ml, 500ml, 1000ml



Product Code	Description
SFH300	Stainless Steel Filter Holder 300ml
SFH500	Stainless Steel Filter Holder 500ml
SFH1000	Stainless Steel Filter Holder 1000ml

SS Filter Holder (Duck Clamp)

- MOC : SS316
- Support screen : SS316
- Stopper : Silicone Rubber
- Clamp : Duck Clamp Anodized Aluminum
- Filter Size : 25mm, 47mm, 90 mm
- Filtration area : 2.1 cm², 9.6 cm², 43 cm²
- Funnel Capacity: 300ml, 500ml.



Product Code	Description
FHDC300	SS 316 Filter Holder 300ml
FHDC500	SS 316 Filter Holder 500ml

Hydrofil Filter Holder

- MOC : SS316
- Support screen : SS316
- Stopper : Silicone Rubber
- Clamp : Spin Lock
- Filter Size : 47mm.
- Filtration area : 9.6 cm²
- Funnel Capacity: 300ml



Product Code	Description
HFH 300	Stainless Steel Filter Holder 300ml

Magnetic Filter Holder

- MOC : Polysulfide PPSU
- Support screen : PES
- Stopper : Silicone Rubber
- Clamp : Magnetic
- Filter Size : 47mm.
- Filtration area : 9.6 cm²
- Funnel Capacity: 300ml, 500ml.



Product Code	Description
MAG300	Magnetic Filter Holder 300ml
MAG500	Magnetic Filter Holder 500ml

PSU Filter Holder

- MOC : Polysulfide
- Support screen : PES
- Stopper : Silicone Rubber
- Clamp : Spin Lock
- Filter Size : 47mm
- Filtration area : 9.6 cm².
- Funnel Capacity: 300ml.



Product Code	Description
PSU300	PSU Filter Holder 300ml

Micro Filtration Funnel

- MOC : Polypropylene
- Filter Size : 47mm
- Filtration area : 9.6 cm²
- Funnel Capacity: 250ml
- ETO Sterilization Packing
- 100 pieces per box



Product Code	Description
MFDF100	Micro Filtration Funnel 100ml
MFDF250	Micro Filtration Funnel 250ml

Disposable M-Funnel

- MOC : Polypropylene
- Filter Size : 47mm (GCN)
- Filtration area : 9.6 cm²
- Funnel Capacity: 100ml
- EO Sterilization Packing
- 200 pieces



Product Cde	Description
DFM100	Micro Filtration Funnel 100ml