



Vial Crimper & De Crimping Tools



Vial Crimper & De Crimper

Product description

KLABS Hand crimpers are designed and manufactured using quality materials, and finishing techniques to provide a consistent and dependable seal while maintaining their durability, and long life. The standard hand crimpers are provided with an adjustable stop, allowing the operator to repeatedly achieve a quality crimped seal. All parts that contact the crimp seal are stainless steel or plated.

The crimper jaws are hardened and polished to provide smooth, consistent operation. Laboratory Grade Hand Crimpers are available with or without cushioned grips which will withstand repeated exposure to cleaning and sterilizing solutions, as well as exposure to some autoclave processes. For convenience, Grips can be removed if necessary. Manufactured using a combination of stainless steel, and hard plated, steel components where both corrosion resistance and increased wear and lubricity are required for optimum function. Our Laboratory Grade Hand Crimpers also have a multi-step, micro surface finish and polishing process for components that require additional wear and lubricity characteristics, such as the crimper jaws. The enhanced surface finish provides smoother operation and allows the material during the crimping process to be more fluid, providing a consistently tight and uniformly crimped seal. ROHS compliant..



- Available : Full Stainless steel 304 & Mild steel Grade Material
- Compatible Aluminum Cap : 1) Normal Vials 2) Flip off type (Before Placing Order Once Confirm)

Vial Decapper and Crimper Pliers Set

Product description

Perfect for Decamping 11mm and 20mm Vials and Bottles (Largely used 20mm Crimper)

- 11mm/20mm Twin Action Vial Decapper - Use it to easily remove aluminum crimp tops from vials and bottles
- 20mm Vial Crimper - Use it to easily place aluminum seals on vials and bottles
- Smooth and simple to operate
- Durable, high-quality construction offering consistent, dependable performance
- Receive both a crimper and decamping plier. Best value for your money.



Aluminum Flip off Vial Cap



Aluminum Tear off Vial Cap



SS Vial Crimper	
Product Code	Size
SSC-1001	11
SSC-1002	13
SSC-1003	20
SSC-1004	32
SSC-1005	33
SS Vial Crimper (Flip off)	
Product Code	Size
SSCFO-1001	11
SSCFO-1002	13
SSCFO-1003	20
SSCFO-1004	32
SSCFO-1005	33

MS Vial Crimper	
Product Code	Size
MSC-1001	11
MSC-1002	13
MSC-1003	20
MSC-1004	32
MSC-1005	33
MS Vial Crimper (Flip off)	
Product Code	Size
MSCFO-1001	11
MSCFO-1002	13
MSCFO-1003	20
MSCFO-1004	32
MSCFO-1005	33

SS plier Vial de-crimper	
Product Code	Size
SSPD-1001	20

SS Vial de-crimper	
Product Code	Size
SSDC-1011	11
SSDC-1012	13
SSDC-1013	20
SSDC-1014	32

MS Vial de-crimper	
Product Code	Size
MSDC-1011	11
MSDC-1012	13
MSDC-1013	20
MSDC-1014	32