



## HPLC / GC CONSUMABLES




### Agilent Purge Valve with PTFE Frits

- purge valve, short, with ptfe frit, for 1120, 1220,
- G7110, G7111 pumps, 600 bar max - purge valves



### Bottle Head Assembly Long ( G7120-60007) Including Contains

- GL45 CAP with tube head 3 Port
- 1/8 HPLC tubing 2 mtr
- Peek Fittings 1/8" - Pkt/3
- Last Drop Filter 1/8x 10 um ss 316 Pkt /1
- Tube No. Pkt Pkt/10



### Aluminium G3188-27501

- Ferrule, flexi, inert
- 0.25mm
- Pkt/ 10/pk



### GC Capillary Column Nut 5181-8830

- Ss 316
- 0.25mm
- Capillary Column Nut (fits Agilent injectors
- non-MS detectors hexagonal wrench tight
- Pkt/ 2/pk



### D2 Deuterium Lamp

- D2-200 Deuterium Lamp for T70, T80 and T92 series



### GC Graphite Ferrules Nut 8010-0304

- Ss 316
- 0.8 mm
- Capillary Ferrule Nut (fits Agilent injectors)
- Pkt/ 10/pk