



## Vacuum Filtration Manifold System



### Vacuum Filtration Manifold Complete Kit **D**



**Vacuum Filtration Manifold Complete Kit**  
(Stainless steel Filter holder External Clamp)

In the Kit:

- 3/6 Stage Manifold
- Filter Funnel SS-316 Inbuilt SS Clamp
- Oil Free Vacuum Pump AS45s
- 47mm x 0.45um Cellulose Nitrate Membrane Filter Paper
- Vacuum Trap with Silicon Tube
- Forceps

Part No	Description
SSCVFCIC-03	Filtration Kit (Inbuilt Clamp)
SSCVFCIC-06	Filtration Kit (Inbuilt Clamp)



**Vacuum Filtration Manifold Complete Kit**  
(Stainless steel Filter holder External Clamp)

In the Kit:

- 3/6 Stage Manifold
- Filter Funnel SS-316 Aluminum Duck Clamp
- Oil Free Vacuum Pump AS75s
- 47mm x 0.45um Cellulose Nitrate Membrane Filter Paper
- Vacuum Trap with Silicon Tube
- Forceps

Part No	Description
SSDVFICIC-03	Filtration Kit (Duck Clamp)
SSDVFICIC-06	Filtration Kit (Duck Clamp)

### Vacuum Filtration Manifold **D**



**Vacuum Filtration Manifold Complete Kit**  
(Glass Funnel)

In the Kit:

- 3/6 Stage Manifold
- Filter Funnel SS-316 with External Duck Clamp
- Oil Free Vacuum Pump AS45s
- 47mm x 0.45um Cellulose Nitrate Membrane Filter Paper
- Vacuum Trap with Silicon Tube
- Forceps

Part No	Description
GFVFCIC-03	Filtration Kit ((Duck Clamp)
GFVFCDC-06	Filtration Kit (Duck Clamp)



**Vacuum Filtration Manifold Complete Kit**  
(Microfilt 47)

In the Kit:

- 3 Stage Manifold
- 250 ml Mfdf250 Funnel With Compatible Adaptor
- Oil Free Vacuum Pump AS75s
- 47mm x 0.45um Cellulose Nitrate Membrane Filter Paper
- Vacuum Trap with Silicon Tube
- Forceps

Part No	Description
MFVFCIC-03	Filtration Kit (MFDF250)
MFVFCDC-06	Filtration Kit (MFDF250)