



Precision Balance

**KPG-10000****SPECIFICATIONS:**

Model	KPG-10000
Make	KLABS
Capacity	10000 g
Resolution	0.1 g
Repeatability	± 0.1 g
Linearity	± 0.2 g
Calibration	External 145 x
Pan Size	175 mm LCD
Display	
Response Time	3-4 secs
Units	g, ct, oz
Battery Backup	In-built
Dimensions	205 x 272 x 163 mm
Weight	3 kg

ACCESSORIES: • Adapter