



## Precision Balance (Professional)



**KPG-5001**

### SALIENT FEATURES:

- Auto Calibration with External Weight
- Easy to read LCD Display with backlight
- Multifunction Weighing Units
- Auto Zero tracking, Low Battery Indicator
- Overload Protection Design



**KPG-3001**

- S.S. Weighing Pan, Level indicator, Adjustable Fit
- High Resolution, Quick Weighing, Accurate Result
- In built rechargeable battery for continues use
- With RS-232 port for data output
- Sturdy Aluminum alloy base

### SPECIFICATIONS:

Model	KPG-3001	KPG-5001
Make	ABS	KLABS
Capacity	3000 g	5000 g
Resolution	0.01 g	0.01 g
Repeatability	± 0.01g	± 0.01g
Linearity	± 0.03g	± 0.03g
Minimum Load	0.05 g	0.05 g
Calibration	External	External
Pan Size	Ø 130 mm	170 x 170 mm
Display	LCD Display	
Response Time	3 ~ 4 secs	
Units	kg, g, ct, lb, oz, ozt, dwt, gsm	
Battery Backup	Yes	
Dimensions & Weight	280 x 190 x 80mm & 3 kg	

### ACCESSORIES:

- Standard Weight
- Adapter
- Instruction Manual