



Magnetic Stirrer Hotplate



MODEL: KHMS-2000

SALIENT FEATURES:

- Ultra-Quiet, Brushless Motor
- Clearly view temperature, speed, and time settings at a glance with a bright, easy-to-read digital display.
- Ceramic plate resists corrosion and scratches for long-term durability.
- Premium Ceramic Plate Surface distributes uniform heat with excellent chemical resistance
- Impact-resistant ABS cabinet provides long-lasting structural integrity and professional aesthetics
- Never lose your last settings—automatic memory saves your previous data for faster repeat experiments.

SPECIFICATIONS:

Model	KHMS-2000
Make	KLABS
Display	Digital
Speed (RPM)	100 to 2000 RPM
Resolution	100 RPM
Plate Size (mm)	168 x 168 mm
Plate Type	Ceramic
Max Volume	5000 ml
Max Temperature	RT+5 to 300°C
Temp. Resolution	1°C
Timing Range	0 to 99 hr, 99 mins
Temp. Sensor	PT 100
Temp. Controller	PID
Motor	Brushless Motor
Power	Motor: 15W, Heating: 500 W
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	280 × 182 × 108 mm & 3 kg

STANDARD

ACCESSORIES:

- PT 100 Temp. Sensor
- 2 x Magnetic Stir Bars
- Adapter
- Instruction Manual