

CE KLABS KCM-1999 Conductivity Meter



Description

KLABS KCM-1999 Conductivity Meter is designed for reliable conductivity measurements with compact portable construction suitable for laboratory and field applications.

Features

- LCD Display
- Accurate & Stable Performance
- Compact Design
- Easy Field Use
- Easy Operation

Standard Accessories

- Instrument Carry Case
- Conductivity Cell
- Battery
- User Manual



Technical Specifications

Technical	Specification
Model No.	KCM-1999
Range	0–1999 μ S, 0–19.99 μ S
Resolution	1 μ S
Sensor	Cond. Cell (dip type)
Battery	9V DC
Accuracy	\pm 1% FS
Display	3½ digit LCD
Size	180 × 85 × 30 mm
Weight	250 gm (Approx.)