



Digital Rotational Viscometer Touchscreen



MODEL: KGV-60

SALIENT FEATURES:

- Auto Range function
- 8 Speed selection
- Direct Viscosity reading and temperature display
- Large bright back-light LCD Display
- Viscosity value is accurate to 2 decimals = 0.01 mpa.s
- With automatic scanning function, the instrument can automatically recommend a suitable combination of rotor and speed
- Can automatically display the range of viscosity that can be measured by the selected combination
- Displays shear stress & shear rate

SPECIFICATIONS:

Model	KGVT-60
Make	KLABS
Measuring Range mpa.s	20~2,000,000 mpa.s
Measurement Resolution	0.01 mpa.s (Full scale)
Speed	0.3, 0.6, 1.5, 3, 6, 12, 30, 60 RPM
No of Speed	8 speed settings
No of Standard Spindle	4
Display	Touchscreen
Auto Scan Function	Can automatically scan and recommend the preferred combination of rotor and speed
Temperature Monitoring	Optional - Feature available
Measurement Accuracy	±1% F.S
Repeatability	± 0.5% F.S
Output	RS 232C Interface
Power Supply	AC 220V / 50Hz
Dimension	300 × 300 × 450 mm
Weight	9 kgs

STANDARD

ACCESSORIES:

- 4 spindles
- Spindle Guard Leg
- Power Cord
- Stand Set
- Instruction Manual