



KLABS KDOM-19 Portable DO Meter



Description

KLABS KDOM-19 Portable DO Meter is designed for dissolved oxygen measurement with stable performance and easy field operation.

Features

- LCD Display
- Accurate & Stable Performance
- Compact Design
- Easy Field Use
- Easy Operation

Standard Accessories

- Instrument Carry Case
- D.O Probe
- Battery
- User Manual



LCD
DISPLAY



STABLE
FLAME



SOFTWARE
CONTROL SYS.



LOW NOISE
COMPRESSOR



USB
CONNECTIVITY

Technical Specifications

Technical	Specification
Model No.	KDOM-19
Range	0–19.9 ppm
Resolution	0.1 ppm
Sensor	Amperometric D.O Probe
Battery	9V DC
Accuracy	±2% FS
Temp. Compensation	5–55°C
Display	3½ digit LCD
Size	180 × 85 × 30 mm
Weight	250 gm (Approx.)