




Low Viscosity Adapter for Viscometers



MODEL: KLVA-30

SALIENT FEATURES:

- Ideal for low viscosity materials
- 20 to 30 mL sample is needed for test
- Spindle No. included with adapter
- Especially made for materials less than 20 mpa.s viscosity
- Viscosity Measurement Range:
1 to 2000 mpa.s



STANDARD ACCESSORIES:

- Joint Road
- Frame
- Pipe
- Hook Rod
- Spindle No. 0

Temperature Sensor for Viscometers

TMP-2

compatible for GV-60



TMP-3

compatible for GV-100 & GV-200

