



High Speed Micro Centrifuge Machine

**MODEL: MCX-15**

SALIENT FEATURES:

- Dedicated 3D digital wheels to adjust speed and time
- Programmable operation can store up to 10 programs
- Motor-driven lid locking system for high safety and reliability
- Heavy base eliminates vibration and movement when running
- Powerful centrifugal force of 15100 x g and Max. speed of 15000rpm
- High speed accuracy of ± 20 rpm
- Large LED display for easy reading
- Fast acceleration and braking system
- Comes with Aluminum rotor with rotor capacity: 1.5/2.0/0.2/0.5mL x 12

SPECIFICATIONS:

Model	KMCX-15
Make	KLABS
Speed Range	500~15000 RPM
Speed Increment	100 RPM
Max. RCF	15100, increment: 100
Speed Accuracy	± 20
Rotor Capacity	0.2ml/0.5ml/1.5ml/ 2ml x 12
Timer	Timer: 30s~99 min
Motor	Brushless DC Motor
Safety Features	Door interlock, Over-speed detection,
Noise	Over-temperature detection, Automatic internal diagnosis
Acceleration/Brake Time	Acceleration: 30s, Braking:10s
Power Supply	AC 220V / 50Hz
Dimension & Weight	256 x 225 x 155 mm & 6 kg