



Smart Mini Centrifuge

With timer & variable speed display



MODEL: KMCX-6V

SALIENT FEATURES:

- Suitable for 2.0ml, 1.5ml, 0.2ml centrifuge tubes & 0.2mlx8 PCR tube strips.
- Compact footprint to save valuable space
- Microprocessor controlled for efficient performance
- Variable timer setting from 1 min to 99 min & infinite mode
- 8 place micro centrifuge rotor

SPECIFICATIONS:

Model	KMCX-6V
Make	KLABS
Speed	1000 to 6000 RPM
Max. Load Capacity	1915 x g (Max RCF)
Speed Accuracy	± 100 RPM
Display	Digital
Motor Power	9W
Time	1 to 99 min & infinite
Rotor Capacity	1. 8 x 1.5/2.0 ml
	2. 10 x 0.2/0.5 ml
	3. 4 place PCR strip
Quick Spin	LID Activated ON/OFF
Voltage	24v 1.5amp, AC 220V / 50 Hz
Dimensions & Weight	150 × 130 × 115 mm & 2.5 kg

STANDARD

ACCESSORIES:

- Main Machine
- Rotors
- Instruction Manual