



## Precision Balance



### SALIENT FEATURES:

- Green LED Display
- Front & Back display
- With Rs232 data output
- Removable glass windshield
- Low battery indicator
- In-nuilt rechargeable battery for continuous use
- Piece counting & percentage function
- Multiple weighing units
- Auto zero tracking

### SPECIFICATIONS:

Model	KPG-620	KPG-1050
Make	KLABS	KLABS
Capacity	620 g	1050 g
Resolution	0.01 g	
Repeatability	± 0.01	
Linearity	± 0.02	
Calibration	External Calibration	
Pan Size	Ø 120 mm	
Tare Range	Full	
Operating Temp.	10 to 40°C	
Display	Green LED (Front & back display)	
Response Time	3-4 secs	
Units	g, ct, oz, dwt	
Battery Backup	In-built	
Dimensions	190 x 240 x 100	
Weight	3 kg	

USE  
LABORATORY & INDUSTRIAL  
SOLUTIONS

### ACCESSORIES:

- Windshield Set
- Pan
- Instruction Manual
- Adapter