

CE Rotary Evaporator (with Touchscreen)



MODEL: KRE-2T

SALIENT FEATURES:

- Safety features: Motor over-current protection, earth-leakage protection lifting overload protection
- 5°C above Set temp. Boil-dry protection mode activates
- In-built vacuum controller
- PTFE sealing rings
- Motorized lifting and rising of the condenser clamp

Technical Specifications

Model	KRE-2T
Make	KLABS
Motor Type	Brushless DC Motor
Speed Range	20 to 300 RPM
Motor Power	40W
Display	5" Touchscreen LCD
Clockwise and Anti-clockwise	Yes
Rotary Bottle Volume	1L / \varnothing 131mm opening
Collection Bottle Volume	1L / \varnothing 131mm opening
Heating Temp. Range	RT to 199°C
Heating Control Accuracy	$\pm 1^\circ\text{C}$
Heating Power	1300 W
Condensation Area	0.15m ²
Evaporation Capacity	as
Timer	Yes
Time Setting Range	0-9999 mins
Vacuum Controller	In-built
Power Supply	AC 220v, 50Hz
Dimensions & Weight	700 x 410 x 880 mm, 43 kgs

STANDARD

ACCESSORIES:

- Rotary Bottle
- Collection Bottle
- Condensor
- Adapter
- Instruction Manual