



pH Meter + Stirrer Combo



SALIENT FEATURES:

- Fully Automatic Microprocessor control.
- Ceramic coated work plate
- Adjustable stirring times & speeds
- Automatic temperature compensation ensures accurate readings over the entire range
- Classical benchtop pH meter equipped with large back-light LCD display



MODEL: KPHS-500C

SPECIFICATIONS:

	Model	KPHS-500C
	Make	KLABS
pH	Range	0.00 ~ 14.00 pH
	Accuracy	± 0.01 ph
	Calibration Point	2 points (pH 4.01/7.00/10.01)
mV/ORP	Range	± 1999 mV
	Accuracy	±1mV
Temp.	Range	0~100°C
	Resolution & Accuracy	Resolution: 0.1°C & Accuracy: ±0.5°C
Stirring	Stirring Volume	0~2000 ml
	Stirring Speed	0~1250 RPM
	Top Plate Size & Material	Ceramic
General	Timer Range	0 to 999 minutes
	Temperature Compensation	0 ~ 100°C
	Connector	BNC
	Power Supply	AC 220V / 50 Hz
	Dimensions & Weight	210 x 205 x 75 mm / 1.5 kgs

STANDARD

ACCESSORIES:

- pH Electrode
- Electrode Arm
- pH Buffer Pouches
- Stir Bars
- Instruction Manual