



KLABS KTDS-199 Portable TDS Meter



Description

KLABS KTDS-199 Portable TDS Meter provides accurate total dissolved solids measurement with portable and user-friendly operation for field and laboratory use.

Features

- LCD Display
- Accurate & Stable Performance
- Compact Design
- Easy Field Use
- Easy Operation

Standard Accessories

- Instrument Carry Case
- TDS Cell
- Battery
- User Manual



Technical Specifications

Technical	Specification
Model No.	KTDS-199
Range	0–1999 ppm, 0–19.99 ppt
Resolution	1 ppm
Sensor	TDS Cell (dip type)
Battery	9V DC
Accuracy	±1% FS
Display	3½ digit LCD
Size	180 × 85 × 30 mm
Weight	250 gm (Approx.)