



Ultrasonic Cleaners (Industrial Size)



SALIENT FEATURES:

- Inner Tank & Outer Body is made of Stainless Steel
- Digital Red LED display
- Temperature & Timer function
- Sports ceramic heaters which are efficient
- Drainage & cooling fan
- High performance industrial transducers
- Easy to move with wheels

SPECIFICATIONS:

Model	KUC-1053	KUC-1103
Make	KLABS	KLABS
Capacity	52 L	103 L
Tank Size	620 x 300 x 280 mm	650 x 530 x 300 mm
Ultrasonic Power	1200 W	2000 W
Heating Power	2000 W	4000 W
Tank Draining Valve	Yes	Yes
Wheels (for moving the machine)	Yes	Yes
Timing Setting	Up to 99 mins	Up to 99 mins
Heating Temp.	R.T +5 ~ 100°C	R.T +5 ~ 100°C
Display	LED	LED
Frequency	40 kHz	40 kHz
Lid Material	SS 304	SS 304
Inner Tank Material	SS 304	SS 304
Power Supply	AC 220V / 50 Hz, 3 phase, 4 wire	
Weight	90 kgs	110 kgs

STANDARD

- Main instrument
- Basket

ACCESSORIES:

- Lid
- Instruction Manual