



Smart Vortex Mixer With timer & speed display



MODEL:KVM-4200

SALIENT FEATURES:

- Suitable for 0.2ml to 50ml micro tubes. Also suitable for containers less than 108mm diameter.
- Continuous & touch operation
- Brush-less DC motor for maintenance free long life
- Built in counter balance for steady & safe operation
- 1 to 999 mins time setting
- Programmable Pulse mode

SPECIFICATIONS:

Model	KVM-4200
Make	KLABS
Rotating Speed	300 to 4200 RPM
Shaking Movement	Orbital
Orbital Diameter	4 mm
Motor Power	20 W
Run Type	Continuous / Touch / Speed
Time	1 to 999 minutes & continuous
Increment/Decrement	10 RPM
Pulse Mode	Adjustable
Max. Load Capacity	500 g
Display	Speed & Time
Ambient Temp. Range	5 to 40°C
Ambient Humidity	80%
Voltage	24v 1.5amp, AC 220V / 50 Hz
Dimensions & Weight	150 × 130 × 115 mm & 2.5 kg

STANDARD

ACCESSORIES:

- Main Machine
- Simple test tube head
- Instruction Manual