



Vortex Mixer



MODEL: KVM-28



MODEL: KMINI-30

SPECIFICATIONS:

Model	KVM-28	KMINI-30
Make	KLABS	KLABS
Speed Range	0-2800rpm	0-3000rpm
Shaking movement	Orbital	Orbital
Orbital diameter	4 mm	4 mm
Motor Power	58W	12W
Run type	Continuos / Touch / Speed	Continuos / Touch / Speed
Ambient Temp. Range	5°C to 40°C	5°C to 40°C
Display	None	None
Power	60W	60W
Power Supply	AC 220V / 50Hz	AC 220V / 50Hz
Dimension LxWxH (mm)	240 x 190 x 240 mm	110 x 90 x 80 mm
Weight	2.8kg	0.5kg

STANDARD

ACCESSORIES:

- Standard Head, Flat Rubber Head, Foam Head
- Power Adapter
- Instruction Manual

- Standard Fixed Head
- Power Adapter
- Instruction Manual