


Western Blotting Apparatus (Wet Transfer/Immunoblotting)



MODEL: KWMINI-TRANS

SALIENT FEATURES:

- Efficient and reproducible protein transfer
- Cooling pack absorbs transfer heat
- Runs two blots simultaneously

SPECIFICATIONS:

Model	KWMINI-TRANS
Make	KLABS
Safety Lid	Lid with electrodes & cables
Transfer Cell	Equipped with platinum electrodes
Gel Holder Casette	110 x 90 mm
Foam Pads	Wrap the gel and membrane (5 pieces)
Cooling Units	2 pieces
Blot Transfer Size	100 x 75 mm
Transfer Capacity	1 - 2 pieces
Dimensions & Weight	180 x 130 x 160 mm & 3 kgs

STANDARD

ACCESSORIES:

- Instruction manual
- Sponge liner
- Blot Tank
- 2 Cooling unit
- Blot frame with electrodes
- Blot connector
- Blot Clips
- Blot Lid with electrodes & cable