



## KLABS Digital D.O Meter



### Product Description

KLABS KDDM-27 Digital D.O Meter is a compact and accurate laboratory instrument designed for dissolved oxygen and temperature measurement applications. The instrument features bright LED display, stable performance and easy operation for laboratory and industrial use.

### Features

- D.O and Temperature indication (LT-27)
- Bright LED display
- Compact design
- Accurate and stable result
- Stable Design
- Easy Operation
- Accurate Reading



ACCURATE  
& STABLE



ACCURATE  
TEMPERATURE



STABLE  
DESIGN



EASY  
OPERATION

## Technical Specifications

Technical	Specification
Product Name	KLABS Digital D.O Meter
Model Number	KLABS KDDM-27
Range	0 to 20.0 ppm
Temperature	0 to 100
Resolution	0.1 ppm
Accuracy	0.1ppm $\pm$ 1 count
Temp. Compensation	0°C to 50°C manual
Compensation	$\pm$ 3%
D.O Sensor	Amperometric gold/silver membrane type
Display	3½ digit LED display
Power	230V AC 50 Hz
Weight	2 Kg

## Accessories

Accessories	Details
D.O Probe	Gold/Silver D.O Probe
D.O Membrane	Included
Temperature probe	Included (LT-27)
Rotational stand	Included
Instruction manual	Included
Dust cover	Included