



Ice Flaker Machine



SALIENT FEATURES:

- Fully Automatic Microprocessor control.
- Stainless Steel construction.
- Continuous ice flakes output
- CFC free compressor
- Low water level indication
- Overload protection
- Noiseless operation
- Compact Design and low maintenance

SPECIFICATIONS:

Model	KIFM-30	KIFM-50	KIFM-100	KIFM-150
Make	KLABS	KLABS	KLABS	KLABS
Capacity of Ice Flaker	30 kg / 24 hr	50 kg / 24 hr	100 kg / 24 hr	150 kg / 24 hr
Storage Capacity	10 kg	15 kg	40 kg	40 kg
Storage Bin	ABS Plastic			
Cooling Method	Air (inside R290A refrigerant)			
Ambient Temperature	Below 25°C for optimum working			
Construction	Stainless Steel			
Noise Level	≤52 dB			
Water Supply	3/4" BSP			
Refrigerant	R290A			
Power Supply	AC 220V / 50 Hz			
Input Power	300 W	380 W	520 W	680 W
Dimensions (mm)	330 x 440 x 605	400 x 540 x 700	500 x 611 x 950	500 x 611 x 950
Gross Weight	26 kg	36 kg	58 kg	60 kg