



KLABS Fume Hood With FRP Lining Sheet



Product Description

With remote operated service line. Metallic CRC sheet duly Powder Coated. Interior chemical resistant, durable, FRP or polypropylene integral work walls with Baffle arrangement. Worktop Chemical resistant splash & spillage proof “PP” or FRP or granite. Auxiliary type fume hood with all remote controlled service lines i.e. water supply, two ports for gas/vacuum & LPG gas. Two electrical sockets.

Features

- Metallic CRC sheet body
- Powder coated construction
- Chemical resistant FRP/PP interior
- Remote controlled service lines
- Splash and spillage proof worktop
- Two electrical sockets provided
- Gas/Vacuum/LPG connections available

Optional

- Fume duct PVC per running ft.
- Grid at the back wall of working chamber to hold the glass apparatus
- IQ/OQ/PQ/DQ certificates



Technical Specifications

Parameter	Specification
Product Name	Fume Hood With FRP Lining Sheet
Model Number	KFH-115
Construction	Metallic CRC Sheet
Finish	Powder Coated
Interior	FRP / Polypropylene
Worktop	PP / FRP / Granite
Service Lines	Remote Controlled
Electrical Sockets	Two
Working Table Size (a)	90cm x 75cm x 75cm (3ft x 2½ft x 2½ft)
Working Table Size (b)	120cm x 75cm x 75cm (4ft x 2½ft x 2½ft)
Working Table Size (c)	150cm x 75cm x 75cm (5ft x 2½ft x 2½ft)
Working Table Size (d)	180cm x 75cm x 75cm (6ft x 2½ft x 2½ft)