



KLABS Analog Vortex Mixer



Product Description

The KLABS KAVM-402 Analog Vortex Mixer provides reliable and uniform sample mixing with variable speed regulation. Its compact design and simple operation make it suitable for everyday laboratory applications.

Features

- Variable Speed Control
- Speed Regulation
- Brushless DC Motor
- Compact Design
- Stable Performance
- 4 mm Orbital Motion
- Easy Operation

Application/Uses

- Research Laboratories
- Diagnostic Laboratories
- Chemical Laboratories
- Pharmaceutical Laboratories
- Educational Institutes



LED DISPLAY




STABLE DESIGN



Stable Anti Slip Foot Pad



DC Brushless motor



Touch / Continuous Mode

Technical Specifications

Parameter	Specification
Product Name	Vortex Mixer
Model Number	KAVM-402
Rotating Speed RPM	300 to 4200
Operation Mode	Speed Regulation
Orbital Diameter	4 mm
Max. Load Capacity	500 g
Protection Class	IP21
Rated Voltage	24V, 1.5 Amp
Rated Power	20W
Frequency	50/60Hz
Ambient Temperature	5° to 40°C
Maximum Humidity	80%
Net Weight	500 g (Approx.)