



KLABS KBRM-04 Blood Roller Mixer



Product Description

KLABS KBRM-04 Blood Roller Mixer is designed for gentle and uniform mixing of blood samples, biological fluids and laboratory specimens. The instrument combines rocking and rolling motion to ensure complete sample mixing while maintaining sample integrity. Its compact design and quiet operation make it ideal for continuous laboratory use.

Features

- Rocking and Rolling Action for Complete Mixing
- Compact Four Roller Design
- Continuous Quiet Operation
- Fixed Speed Operation
- Accurate and Stable Performance
- Easy Operation
- Low Noise Level
- Suitable for 24-Hour Continuous Working

Application/Uses

- Blood Sample Mixing
- Clinical Laboratories
- Hematology Laboratories
- Biological Sample Preparation
- Research Laboratories
- Diagnostic Centers
- Sample Homogenization Applications



FIXED SPEED
25-30 RPM



ACCURATE
& STABLE



EASY
OPERATION



LOW NOISE



24 HOURS
WORKING

Technical Specifications

Parameter	Details
Product Name	Blood Roller Mixer
Model Number	KBRM-04
Roller Size	230 mm
Controls	Analog
Number of Rollers	4
Speed	33 RPM (Approx.)
Maximum Load	5 Test Tubes
Timer	Not Available
Power Supply	230V $\pm$ 10% AC, 50 Hz

Standard Accessories

- Power Cord
- Instruction Manual
- Roller Platform Assembly

Advantages

- Gentle Mixing Without Sample Damage
- Uniform Sample Suspension
- Compact Laboratory Footprint
- Maintenance-Free Operation
- Reliable Continuous Performance
- Easy Sample Loading and Removal