



KLABS Mini Vortex Mixer



Product Description



The KLABS KMVM-142 Mini Vortex Mixer is a lightweight and compact mixer designed for quick sample mixing. It provides fixed-speed operation and is ideal for laboratories requiring fast and convenient mixing of small-volume samples.

Features

- Fixed Speed Operation
- Compact & Portable
- One-Touch Mixing
- Stable Design
- Low Maintenance
- Suitable for Routine Laboratory Use

Application/Uses

- Clinical Laboratories
- Biotechnology Laboratories
- Research Laboratories
- Educational Institutes
- Sample Preparation



Technical Specifications

Parameter	Specification
Product Name	Mini Vortex Mixer
Model Number	KMVM-142
Rotating Speed RPM	4200
Orbital Diameter	4 mm
Max. Load Capacity	500 g
Protection Class	IP21
Rated Voltage	24V, 1.5 Amp
Rated Power	20W
Frequency	50/60Hz
Ambient Temperature	5° to 40°C
Maximum Humidity	80%
Net Weight	500 g (Approx.)