



KLABS Portable Haemoglobin Meter



Product Description

KLABS KPHM-25 Portable Haemoglobin Meter is a compact and easy-to-use device designed for quick haemoglobin measurement. It provides accurate and stable results with simple operation, making it suitable for clinical, laboratory and field use.

Features

- Compact body design
- Accurate and stable result
- 1000 test results
- LCD display
- 24 hours working
- Easy operation

Application/Uses

- Quick haemoglobin testing
- Anaemia screening
- Clinic and laboratory use
- Health camp and field testing

Accessories

- HB meter
- Canister of test strips
- Code chip
- Control chip
- Capillary transfer tube/droppe
- Lancing device
- Sterile lancets
- AAA batteries



LCD
DISPLAY



ACCURATE
& STABLE



24 HOURS
WORKING



COMPACT
DESIGN



EASY
OPERATION

Technical Specifications

Specification	Details
Product Name	Portable Haemoglobin Meter
Model Number	KPHM-25
Testing Strips	10 strips
Measuring Range	5-25 gm/dl
Memory	1000 test results
Calibration	Auto calibration