



Mini Centrifuge

With timer & speed display

**MODEL: KMCX-12**

SALIENT FEATURES:

- Brushless motor for long life
- 2 adapters included to use PCR & MCT tubes
- Automatic motor shut down on opening the cover
- Compact design
- LCD display and microprocessor control

SPECIFICATIONS:

Model	KMCX-12
Make	KLABS
Speed Range	100 to 12000 RPM
Speed Adjustment	Digital
Speed Increament	10 RPM
Speed Control Accuracy	±3%
Relative centrifugal force	7245 g
Timer	1 to 99 minutes
Motor Type	Bruchless DC Motor
Motor Power	48 W
Display	LCD
Noise	51-53 dBA
Rotor	MCT 8*2mL, PCR 16*0.2mL
Tubes which can be used	2/1.5/0.5/0.2m8*2mL,16*0.2mL
Adapter	0.5ml, 0.2ml Adapter
Dimensions & Weight	170 x 146 x 119 mm & 1 kg

STANDARD

ACCESSORIES:

- 8 Hole Rotor
- 16 Hole Rotor
- 0.5ml Adapter
- 0.5ml Adapter
- Power Adapter
- Instruction Manual